

ZNF134 Polyclonal Antibody

Catalog No: #41542

Package Size: #41542-1 50ul #41542-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ZNF134 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Mk
Specificity	ZNF134 Polyclonal Antibody detects endogenous levels of ZNF134 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ZNF134.
Target Name	ZNF134
Other Names	ZNF134; Zinc finger protein 134
Accession No.	Swiss-Prot: P52741NCBI Gene ID: 7693
Uniprot	P52741
GeneID	7693;
SDS-PAGE MW	36kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of K562 cells using ZNF134 Polyclonal Antibody

Note: This product is for in vitro research use only