

ZNF695 Polyclonal Antibody

Catalog No: #41544

Package Size: #41544-1 50ul #41544-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZNF695 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	ZNF695 Polyclonal Antibody detects endogenous levels of ZNF695 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ZNF695.
Target Name	ZNF695
Other Names	ZNF695; Zinc finger protein 695; Zinc finger protein SBZF3
Accession No.	Swiss-Prot: Q8IW36NCBI Gene ID: 57116
Uniprot	Q8IW36
GeneID	57116;
SDS-PAGE MW	20kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 293T cells using ZNF695 Polyclonal Antibody

Note: This product is for in vitro research use only