

ENaC γ Polyclonal Antibody

Catalog No: #41547

Package Size: #41547-1 50ul #41547-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ENaC γ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	ENaC γ Polyclonal Antibody detects endogenous levels of ENaC γ protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ENaC γ .
Target Name	ENaC γ
Other Names	SCNN1G; Amiloride-sensitive sodium channel subunit gamma; Epithelial Na(+) channel subunit gamma; ENaCG; Gamma-ENaC; Gamma-NaCH; Nonvoltage-gated sodium channel 1 subunit gamma; SCNEG
Accession No.	Swiss-Prot: P51170NCBI Gene ID: 6340
Uniprot	P51170
GeneID	6340;
SDS-PAGE MW	80kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images

Western Blot analysis of A549 cells using ENaC γ Polyclonal Antibody

Note: This product is for in vitro research use only