

GAPDH Polyclonal Antibody

Catalog No: #41549

Package Size: #41549-1 50ul #41549-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	GAPDH Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA IF
Species Reactivity	Hu Ms Rt
Specificity	GAPDH Polyclonal Antibody detects endogenous levels of GAPDH protein.
Immunogen Description	Synthesized peptide derived from the N-terminal region of human GAPDH.
Target Name	GAPDH
Other Names	GAPDH; GAPD; CDABP0047; OK/SW-cl.12; Glyceraldehyde-3-phosphate dehydrogenase; GAPDH; Peptidyl-cysteine S-nitrosylase GAPDH
Accession No.	Swiss-Prot: P04406NCBI Gene ID: 2597
Uniprot	P04406
GeneID	2597;
SDS-PAGE MW	36kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of NIH-3T3 HepG2 MCF7 cells using GAPDH Polyclonal Antibody

Note: This product is for in vitro research use only