

VEGF-A Polyclonal Antibody

Catalog No: #41552

Package Size: #41552-1 50ul #41552-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	VEGF-A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	VEGF-A Polyclonal Antibody detects endogenous levels of VEGF-A protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human VEGF-A.
Target Name	VEGF-A
Other Names	VEGFA; VEGF; Vascular endothelial growth factor A; VEGF-A; Vascular permeability factor; VPF
Accession No.	Swiss-Prot: P15692NCBI Gene ID: 7422
Uniprot	P15692
GeneID	7422;
SDS-PAGE MW	27kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 22RV-1 SW480 BT474 K562 mouse liver cells using VEGF-A Polyclonal Antibody

Note: This product is for in vitro research use only