

## IL-10 Polyclonal Antibody

Catalog No: #41553

Package Size: #41553-1 50ul #41553-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	IL-10 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA IF
Species Reactivity	Hu
Specificity	IL-10 Polyclonal Antibody detects endogenous levels of IL-10 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human IL-10.
Target Name	IL-10
Other Names	IL10; Interleukin-10; IL-10; Cytokine synthesis inhibitory factor; CSIF
Accession No.	Swiss-Prot: P22301NCBI Gene ID: 3586
Uniprot	P22301
GeneID	3586;
SDS-PAGE MW	20kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

IF 1:100-300


ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of THP-1 Romas HepG2 SW480 SHSY5Y Mouse liver cells using IL-10 Polyclonal Antibody



Immunofluorescence analysis of A549.

---

Note: This product is for in vitro research use only