

## Sg III Polyclonal Antibody

Catalog No: #41557

Package Size: #41557-1 50ul #41557-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Sg III Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Rt
Specificity	Sg III Polyclonal Antibody detects endogenous levels of Sg III protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Sg III.
Target Name	Sg III
Other Names	SCG3; Secretogranin-3; Secretogranin III; SgIII
Accession No.	Swiss-Prot: Q8WXD2NCBI Gene ID: 29106
Uniprot	Q8WXD2
GeneID	29106;
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.

## Images



Western Blot analysis of mouse heart mouse kidney cells using Sg III Polyclonal Antibody

---

Note: This product is for in vitro research use only