

Cathepsin B Rabbit mAb

Catalog No: #49515



Package Size: #49515-1 50ul #49515-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cathepsin B Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	JA11-02
Purification	ProA affinity purified
Applications	WB, ICC/IF, IHC
Species Reactivity	Hu
Immunogen Description	recombinant protein
Other Names	Amyloid precursor protein secretase antibody APP secretase antibody APPS antibody CATB_HUMAN antibody Cathepsin B heavy chain antibody Cathepsin B1 antibody CathepsinB antibody CPSB antibody CTSB antibody cysteine protease antibody OTTHUMP00000116009 antibody OTTHUMP00000229510 antibody OTTHUMP00000229511 antibody OTTHUMP00000229512 antibody OTTHUMP00000229514 antibody OTTHUMP00000229515 antibody OTTHUMP00000229516 antibody Preprocathepsin B antibody
Accession No.	Swiss-Prot#:P07858
Uniprot	P07858
GeneID	1508;
Calculated MW	38 kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

WB: 1:500-1:1,000 IHC: 1:50-1:200 ICC: 1:50-1:200

Background

The cathepsin family of proteolytic enzymes contains several diverse classes of proteases. The cysteine protease class comprises cathepsins B, L, H, K, S, and O. The aspartyl protease class is composed of cathepsins D and E. Cathepsin G is in the serine protease class. Most cathepsins are lysosomal and each is involved in cellular metabolism, participating in various events such as peptide biosynthesis and protein degradation. Cathepsin B is expressed in luminal epithelial cells, indicating that cathepsin B is a marker for secretory cell death.

References

Note: This product is for in vitro research use only