

PKC γ Polyclonal Antibody

Catalog No: #41563

Package Size: #41563-1 50ul #41563-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PKC γ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	PKC γ Polyclonal Antibody detects endogenous levels of PKC γ protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human PKC γ .
Target Name	PKC γ
Other Names	PRKCG; PKCG; Protein kinase C gamma type; PKC-gamma
Accession No.	Swiss-Prot: P05129NCBI Gene ID: 5582
Uniprot	P05129
GeneID	5582;
SDS-PAGE MW	78kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of A549 THP-1 K562 mouse brain mouse skeletal muscle mouse testis cells using PKC γ Polyclonal Antibody

Note: This product is for in vitro research use only