

FAF1 Polyclonal Antibody

Catalog No: #41564

Package Size: #41564-1 50ul #41564-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	FAF1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	FAF1 Polyclonal Antibody detects endogenous levels of FAF1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human FAF1.
Target Name	FAF1
Other Names	FAF1; UBXD12; UBXN3A; CGI-03; FAS-associated factor 1; hFAF1; UBX domain-containing protein 12; UBX domain-containing protein 3A
Accession No.	Swiss-Prot: Q9UNN5NCBI Gene ID: 11124
Uniprot	Q9UNN5
GeneID	11124;
SDS-PAGE MW	73kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of K562 cells using FAF1 Polyclonal Antibody

Note: This product is for in vitro research use only