

BNIP-3 Polyclonal Antibody

Catalog No: #41567

Package Size: #41567-1 50ul #41567-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BNIP-3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	BNIP-3 Polyclonal Antibody detects endogenous levels of BNIP-3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human BNIP-3.
Target Name	BNIP-3
Other Names	BNIP3; NIP3; BCL2/adenovirus E1B 19 kDa protein-interacting protein 3
Accession No.	Swiss-Prot: Q12983NCBI Gene ID: 664
Uniprot	Q12983
GeneID	664;
SDS-PAGE MW	22kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of mouse eye cells using BNIP-3 Polyclonal Antibody

Note: This product is for in vitro research use only