

Catalase Polyclonal Antibody

Catalog No: #41568

Package Size: #41568-1 50ul #41568-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Catalase Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Catalase Polyclonal Antibody detects endogenous levels of Catalase protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Catalase.
Target Name	Catalase
Other Names	CAT; Catalase
Accession No.	Swiss-Prot: P04040NCBI Gene ID: 847
Uniprot	P04040
GeneID	847;
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HT29 A549 cells using Catalase Polyclonal Antibody

Note: This product is for in vitro research use only