

FGF-12 Polyclonal Antibody

Catalog No: #41572

Package Size: #41572-1 50ul #41572-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FGF-12 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	FGF-12 Polyclonal Antibody detects endogenous levels of FGF-12 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human FGF-12.
Target Name	FGF-12
Other Names	FGF12; FGF12B; FHF1; Fibroblast growth factor 12; FGF-12; Fibroblast growth factor homologous factor 1; FHF-1; Myocyte-activating factor
Accession No.	Swiss-Prot: P61328NCBI Gene ID: 2257
Uniprot	P61328
GeneID	2257;
SDS-PAGE MW	27kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of 22RV-1 BT474 SHSY5Y cells using FGF-12 Polyclonal Antibody

Note: This product is for in vitro research use only