

HSP40-4 Polyclonal Antibody

Catalog No: #41573

Package Size: #41573-1 50ul #41573-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	HSP40-4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA IHC-p
Species Reactivity	Hu Ms Rt
Specificity	HSP40-4 Polyclonal Antibody detects endogenous levels of HSP40-4 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human HSP40-4.
Target Name	HSP40-4
Other Names	DNAJA1; DNAJ2; HDJ2; HSJ2; HSPF4; DnaJ homolog subfamily A member 1; DnaJ protein homolog 2; HSDJ; Heat shock 40 kDa protein 4; Heat shock protein J2; HSJ-2; Human DnaJ protein 2; hDj-2
Accession No.	Swiss-Prot: P31689NCBI Gene ID: 3301
Uniprot	P31689
GeneID	3301;
SDS-PAGE MW	44kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

IHC-p: 1:100-300


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of mouse liver mouse kidney mouse brain cells using HSP40-4 Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded rat-testis, antibody was diluted at 1:100

Note: This product is for in vitro research use only