

Trk C Polyclonal Antibody

Catalog No: #41576

Package Size: #41576-1 50ul #41576-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Trk C Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Trk C Polyclonal Antibody detects endogenous levels of Trk C protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Trk C.
Target Name	Trk C
Other Names	NTRK3; TRKC; NT-3 growth factor receptor; GP145-TrkC; Trk-C; Neurotrophic tyrosine kinase receptor type 3; TrkC tyrosine kinase
Accession No.	Swiss-Prot: Q16288NCBI Gene ID: 4916
Uniprot	Q16288
GeneID	4916;
SDS-PAGE MW	95kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using Trk C Polyclonal Antibody

Note: This product is for in vitro research use only