

FGF-6 Polyclonal Antibody

Catalog No: #41580

Package Size: #41580-1 50ul #41580-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FGF-6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	FGF-6 Polyclonal Antibody detects endogenous levels of FGF-6 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human FGF-6.
Target Name	FGF-6
Other Names	FGF6; HST2; HSTF2; Fibroblast growth factor 6; FGF-6; Heparin secretory-transforming protein 2; HST-2; HSTF-2; Heparin-binding growth factor 6; HBGF-6
Accession No.	Swiss-Prot: P10767NCBI Gene ID: 2251
Uniprot	P10767
GeneID	2251;
SDS-PAGE MW	22kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of 293 HepG2 PC-3 M21 cells using FGF-6 Polyclonal Antibody

Note: This product is for in vitro research use only