

Angptl4 Polyclonal Antibody

Catalog No: #41583

Package Size: #41583-1 50ul #41583-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angptl4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Angptl4 Polyclonal Antibody detects endogenous levels of Angptl4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Angptl4.
Target Name	Angptl4
Other Names	ANGPTL4; ARP4; HFARP; PGAR; PP1158; PSEC0166; Angiopoietin-related protein 4; Angiopoietin-like protein 4; Hepatic fibrinogen/angiopoietin-related protein; HFARP
Accession No.	Swiss-Prot: Q9BY76NCBI Gene ID: 51129
Uniprot	Q9BY76
GeneID	51129;
SDS-PAGE MW	45kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using Angptl4 Polyclonal Antibody

Note: This product is for in vitro research use only