

## IL-16 Polyclonal Antibody

Catalog No: #41584

Package Size: #41584-1 50ul #41584-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

|                       |   |
|-----------------------|---|
| Product Name          | IL-16 Polyclonal Antibody   |
| Host Species          | Rabbit  |
| Clonality             | Polyclonal  |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB ELISA  |
| Species Reactivity    | Hu  |
| Specificity           | IL-16 Polyclonal Antibody detects endogenous levels of IL-16 protein.   |
| Immunogen Type        | peptide   |
| Immunogen Description | Synthesized peptide derived from the C-terminal region of human IL-16.  |
| Target Name           | IL-16   |
| Other Names           | IL16; Pro-interleukin-16  |
| Accession No.         | Swiss-Prot: Q14005NCBI Gene ID: 3603  |
| Uniprot               | Q14005  |
| GeneID                | 3603;   |
| SDS-PAGE MW           | 140kd   |
| Concentration         | 1mg/ml  |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| Storage               | Store at -20°C/1 year   |


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of 293 HepG2 BXPC-3 Raji mouse lung mouse brain mouse spleen cells using IL-16 Polyclonal Antibody

---

Note: This product is for in vitro research use only