

FGF-18 Polyclonal Antibody

Catalog No: #41590

Package Size: #41590-1 50ul #41590-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	FGF-18 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	FGF-18 Polyclonal Antibody detects endogenous levels of FGF-18 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human FGF-18.
Target Name	FGF-18
Other Names	FGF18; Fibroblast growth factor 18; FGF-18; zFGF5
Accession No.	Swiss-Prot: O76093NCBI Gene ID: 8817
Uniprot	O76093
GeneID	8817;
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of SW480 SKOV3 BT474 mouse heart cells using FGF-18 Polyclonal Antibody

Note: This product is for in vitro research use only