

Aldolase C Polyclonal Antibody

Catalog No: #41593

Package Size: #41593-1 50ul #41593-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Aldolase C Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Aldolase C Polyclonal Antibody detects endogenous levels of Aldolase C protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Aldolase C.
Target Name	Aldolase C
Other Names	ALDOC; ALDC; Fructose-bisphosphate aldolase C; Brain-type aldolase
Accession No.	Swiss-Prot: P09972NCBI Gene ID: 230
Uniprot	P09972
GeneID	230;
SDS-PAGE MW	39kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of mouse liver mouse spinal cord cells using Aldolase C Polyclonal Antibody

Note: This product is for in vitro research use only