

FGF-5 Polyclonal Antibody

Catalog No: #41599

Package Size: #41599-1 50ul #41599-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	FGF-5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	FGF-5 Polyclonal Antibody detects endogenous levels of FGF-5 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human FGF-5.
Target Name	FGF-5
Other Names	FGF5; Fibroblast growth factor 5; FGF-5; Heparin-binding growth factor 5; HBGF-5; Smag-82
Accession No.	Swiss-Prot: P12034NCBI Gene ID: 2250
Uniprot	P12034
GeneID	2250;
SDS-PAGE MW	30kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of mouse liver cells using FGF-5 Polyclonal Antibody

Note: This product is for in vitro research use only