

PKC gamma Rabbit mAb

Catalog No: #49912

Package Size: #49912-1 50ul #49912-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PKC gamma Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	JG69-32
Purification	ProA affinity purified
Applications	WB,IHC
Species Reactivity	Hu, Ms, Rt
Immunogen Description	Synthetic peptide within Human PKC gamma aa 350-450.
Other Names	KPCG_HUMAN antibody MGC57564 antibody OTTHUMP0000067291 antibody PKC-gamma antibody PKCC antibody PKCG antibody PRKCG antibody Protein kinase C gamma antibody Protein kinase C gamma polypeptide antibody Protein kinase C gamma type antibody Protein kinase C, gamma antibody SCA 14 antibody SCA14 antibody
Accession No.	Swiss-Prot#:P05129
Uniprot	P05129
GeneID	5582;
Calculated MW	78 kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C


Application Details

WB: 1:500-1:2,000 IHC: 1:50-1:200


Images



Western blot analysis of PKC gamma on A549 cell lysates using anti-PKC gamma antibody at 1/1,000 dilution.



Immunohistochemical analysis of paraffin-embedded rat cerebellum tissue using anti-PKC gamma antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded mouse cerebellum tissue using anti-PKC gamma antibody. Counter stained with hematoxylin.

References

Note: This product is for in vitro research use only