

ASIC3 Polyclonal Antibody

Catalog No: #41609

Package Size: #41609-1 50ul #41609-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ASIC3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	ASIC3 Polyclonal Antibody detects endogenous levels of ASIC3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ASIC3.
Target Name	ASIC3
Other Names	ASIC3; ACCN3; SLNAC1; TNAC1; Acid-sensing ion channel 3; ASIC3; hASIC3; Amiloride-sensitive cation channel 3; Neuronal amiloride-sensitive cation channel 3; Testis sodium channel 1; hTNaC1
Accession No.	Swiss-Prot: Q9UHC3NCBI Gene ID: 9311
Uniprot	Q9UHC3
GeneID	9311;
SDS-PAGE MW	58kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using ASIC3 Polyclonal Antibody

Note: This product is for in vitro research use only