

IL-1 $\beta$  Polyclonal Antibody

Catalog No: #41610

Package Size: #41610-1 50ul #41610-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	IL-1 $\beta$ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA IHC-p IF/ICC
Species Reactivity	Hu Ms Rt
Specificity	IL-1 $\beta$ Polyclonal Antibody detects endogenous levels of IL-1 $\beta$ protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human IL-1 $\beta$ .
Target Name	IL-1 $\beta$
Other Names	IL1B; IL1F2; Interleukin-1 beta; IL-1 beta; Catabolin
Accession No.	Swiss-Prot: P01584NCBI Gene ID: 3553
Uniprot	P01584
GeneID	3553;
SDS-PAGE MW	17kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.


IF 1:100-300

IHC-p: 1:100-300

ELISA: 1/20000.

Not yet tested in other applications.

## Images

Western Blot analysis of ES-2 U87 PC-3 SHSY5Y mouse liver cells using IL-1 $\beta$  Polyclonal Antibody

---

Note: This product is for in vitro research use only