

ADM Polyclonal Antibody

Catalog No: #41611

Package Size: #41611-1 50ul #41611-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADM Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Rt
Specificity	ADM Polyclonal Antibody detects endogenous levels of ADM protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ADM.
Target Name	ADM
Other Names	ADM; AM; ADM
Accession No.	Swiss-Prot: P35318NCBI Gene ID: 133
Uniprot	P35318
GeneID	133;
SDS-PAGE MW	20kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of mouse kidney mouse lung mouse heart cells using ADM Polyclonal Antibody

Note: This product is for in vitro research use only