

Annexin II Polyclonal Antibody

Catalog No: #41612

Package Size: #41612-1 50ul #41612-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Annexin II Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Annexin II Polyclonal Antibody detects endogenous levels of Annexin II protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Annexin II.
Target Name	Annexin II
Other Names	ANXA2; ANX2; ANX2L4; CAL1H; LPC2D; Annexin A2; Annexin II; Annexin-2; Calpactin I heavy chain; Calpactin-1 heavy chain; Chromobindin-8; Lipocortin II; Placental anticoagulant protein IV; PAP-IV; Protein I; p36
Accession No.	Swiss-Prot: P07355NCBI Gene ID: 302
Uniprot	P07355
GeneID	302;
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa HepG2 cells using Annexin II Polyclonal Antibody

Note: This product is for in vitro research use only