

Lipocalin-1 Polyclonal Antibody

Catalog No: #41627

Package Size: #41627-1 50ul #41627-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lipocalin-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Lipocalin-1 Polyclonal Antibody detects endogenous levels of Lipocalin-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Lipocalin-1.
Target Name	Lipocalin-1
Other Names	LCN1; VEGP; Lipocalin-1; Tear lipocalin; Tlc; Tear prealbumin; TP; Von Ebner gland protein; VEG protein
Accession No.	Swiss-Prot: P31025NCBI Gene ID: 3933
Uniprot	P31025
GeneID	3933;
SDS-PAGE MW	21kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of K562 cells using Lipocalin-1 Polyclonal Antibody

Note: This product is for in vitro research use only