

# mHTF culture medium

Catalog No: #MSE-3018



Package Size: #MSE-3018-1 10mL/bottle #MSE-3018-2 30mL/bottle

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	mHTF culture medium
Storage	The kit is stored at 2~8°C in the dark and can be stored for 3 months. It is valid for 1 week after opening. Avoid repeated freezing and thawing.

## Product Description

Composition: Purpose: This product uses HEPES buffer system to maintain the stability of environment (PH value), and is for adjustable medium of all gametes and pre-implantation embryo operation. Quality Index: 1. pH: 7.2 0.4

2. Osmotic pressure: 30020 mOsm/kg

3. Mouse embryo test MEA: Blastocyst formation rate 80% (1-cell was cultured for 96 hours.)

4. Sterilization: filtration method (SAL 10<sup>-6</sup>) Sample Requirements: Using Protocol: The culture medium can only be used in the air environment. There is no need to balance before use. Attentions: This product is for research use only.

Note: This product is for in vitro research use only