

Layilin Polyclonal Antibody

Catalog No: #41630

Package Size: #41630-1 50ul #41630-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Layilin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Layilin Polyclonal Antibody detects endogenous levels of Layilin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Layilin.
Target Name	Layilin
Other Names	LAYN; Layilin
Accession No.	Swiss-Prot: Q6UX15NCBI Gene ID: 143903
Uniprot	Q6UX15
GeneID	143903;
SDS-PAGE MW	43kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of HuvEc K562 cells using Layilin Polyclonal Antibody

Note: This product is for in vitro research use only