

DPYD Polyclonal Antibody

Catalog No: #41635

Package Size: #41635-1 50ul #41635-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DPYD Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	DPYD Polyclonal Antibody detects endogenous levels of DPYD protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human DPYD.
Target Name	DPYD
Other Names	DPYD; Dihydropyrimidine dehydrogenase [NADP(+)]; DHPDHase; DPD; Dihydrothymine dehydrogenase; Dihydrouracil dehydrogenase
Accession No.	Swiss-Prot: Q12882NCBI Gene ID: 1806
Uniprot	Q12882
GeneID	1806;
SDS-PAGE MW	120kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of HepG2 cells using DPYD Polyclonal Antibody

Note: This product is for in vitro research use only