

VHR Polyclonal Antibody

Catalog No: #41636

Package Size: #41636-1 50ul #41636-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	VHR Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	VHR Polyclonal Antibody detects endogenous levels of VHR protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human VHR.
Target Name	VHR
Other Names	DUSP3; VHR; Dual specificity protein phosphatase 3; Dual specificity protein phosphatase VHR; Vaccinia H1-related phosphatase; VHR
Accession No.	Swiss-Prot: P51452NCBI Gene ID: 1845
Uniprot	P51452
GeneID	1845;
SDS-PAGE MW	21kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using VHR Polyclonal Antibody

Note: This product is for in vitro research use only