

Neuropilin-2 Polyclonal Antibody

Catalog No: #41638

Package Size: #41638-1 50ul #41638-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neuropilin-2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Neuropilin-2 Polyclonal Antibody detects endogenous levels of Neuropilin-2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Neuropilin-2.
Target Name	Neuropilin-2
Other Names	NRP2; VEGF165R2; Neuropilin-2; Vascular endothelial cell growth factor 165 receptor 2
Accession No.	Swiss-Prot: O60462NCBI Gene ID: 8828
Uniprot	O60462
GeneID	8828;
SDS-PAGE MW	100kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using Neuropilin-2 Polyclonal Antibody

Note: This product is for in vitro research use only