

SPTLC1 Polyclonal Antibody

Catalog No: #41640

Package Size: #41640-1 50ul #41640-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SPTLC1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	SPTLC1 Polyclonal Antibody detects endogenous levels of SPTLC1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human SPTLC1.
Target Name	SPTLC1
Other Names	SPTLC1; LCB1; Serine palmitoyltransferase 1; Long chain base biosynthesis protein 1; LCB 1; Serine-palmitoyl-CoA transferase 1; SPT 1; SPT1
Accession No.	Swiss-Prot: O15269NCBI Gene ID: 10558
Uniprot	O15269
GeneID	10558;
SDS-PAGE MW	52kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of HepG2 cells using SPTLC1 Polyclonal Antibody

Note: This product is for in vitro research use only