

## LASP-1 Polyclonal Antibody

Catalog No: #41642

Package Size: #41642-1 50ul #41642-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	LASP-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	LASP-1 Polyclonal Antibody detects endogenous levels of LASP-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human LASP-1.
Target Name	LASP-1
Other Names	LASP1; MLN50; LIM and SH3 domain protein 1; LASP-1; Metastatic lymph node gene 50 protein; MLN 50
Accession No.	Swiss-Prot: Q14847NCBI Gene ID: 3927
Uniprot	Q14847
GeneID	3927;
SDS-PAGE MW	28kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.

## Images



Western Blot analysis of HepG2 cells using LASP-1 Polyclonal Antibody

---

Note: This product is for in vitro research use only