

PEPCK-C Polyclonal Antibody

Catalog No: #41644

Package Size: #41644-1 50ul #41644-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PEPCK-C Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PEPCK-C Polyclonal Antibody detects endogenous levels of PEPCK-C protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human PEPCK-C.
Target Name	PEPCK-C
Other Names	PCK1; PEPCK1; Phosphoenolpyruvate carboxykinase, cytosolic [GTP]; PEPCK-C; Phosphoenolpyruvate carboxylase
Accession No.	Swiss-Prot: P35558NCBI Gene ID: 5105
Uniprot	P35558
GeneID	5105;
SDS-PAGE MW	65kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of HeLa MCF7 cells using PEPCK-C Polyclonal Antibody

Note: This product is for in vitro research use only