

CD1C Polyclonal Antibody

Catalog No: #41646

Package Size: #41646-1 50ul #41646-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD1C Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CD1C Polyclonal Antibody detects endogenous levels of CD1C protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD1C.
Target Name	CD1C
Other Names	CD1C; T-cell surface glycoprotein CD1c; CD1c
Accession No.	Swiss-Prot: P29017NCBI Gene ID: 911
Uniprot	P29017
GeneID	911;
SDS-PAGE MW	37kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of K562 cells using CD1C Polyclonal Antibody

Note: This product is for in vitro research use only