

## CD1E Polyclonal Antibody

Catalog No: #41647

Package Size: #41647-1 50ul #41647-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	CD1E Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CD1E Polyclonal Antibody detects endogenous levels of CD1E protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human CD1E.
Target Name	CD1E
Other Names	CD1E; T-cell surface glycoprotein CD1e, membrane-associated; hCD1e; R2G1; CD1e
Accession No.	Swiss-Prot: P15812NCBI Gene ID: 913
Uniprot	P15812
GeneID	913;
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of K562 cells using CD1E Polyclonal Antibody

---

Note: This product is for in vitro research use only