

IL-2R α Polyclonal Antibody

Catalog No: #41650

Package Size: #41650-1 50ul #41650-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | IL-2R α Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications | WB ELISA |
| Species Reactivity | Hu |
| Specificity | IL-2R α Polyclonal Antibody detects endogenous levels of IL-2R α protein. |
| Immunogen Type | peptide |
| Immunogen Description | Synthesized peptide derived from the Internal region of human IL-2R α . |
| Target Name | IL-2R α |
| Other Names | IL2RA; Interleukin-2 receptor subunit alpha; IL-2 receptor subunit alpha; IL-2-RA; IL-2R subunit alpha; IL2-RA; TAC antigen; p55; CD25 |
| Accession No. | Swiss-Prot: P01589NCBI Gene ID: 3559 |
| Uniprot | P01589 |
| GeneID | 3559; |
| SDS-PAGE MW | 30kd |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images

Western Blot analysis of K562 cells using IL-2R α Polyclonal Antibody

Note: This product is for in vitro research use only