

Mouse Tumor necrosis factor-like ligand 1 (TL1) ELISA Kit

Catalog No: #EK6200

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #EK6200-1 48T #EK6200-2 96T

Description

Product Name	Mouse Tumor necrosis factor-like ligand 1 (TL1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse (<i>Mus musculus</i>)
Other Names	MGC129934; MGC129935; TL1; TL1A; VEGI; VEGI192A; TNF ligand-related molecule 1 TNF superfamily ligand TL1A vascular endothelial cell growth inhibitor vascular endothelial growth inhibitor-192A
Accession No.	Q5UBV8
Uniprot	Q5UBV8
GeneID	326623;
Storage	<p>The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.</p> <p>The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days, and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).</p>

Application Details

Detect Range:62.5-4000 pg/mL

Sensitivity:26.8 pg/mL

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:SandwichTest principle:This assay employs a two-site sandwich ELISA to quantitate TNFSF15 in samples. An antibody specific for TNFSF15 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyTNFSF15 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TNFSF15 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TNFSF15 bound in the initial step. The color development is stopped and the intensity of the color is measured.

Note: This product is for in vitro research use only