

CXCR-3 Polyclonal Antibody

Catalog No: #41683

Package Size: #41683-1 50ul #41683-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CXCR-3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CXCR-3 Polyclonal Antibody detects endogenous levels of CXCR-3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CXCR-3.
Target Name	CXCR-3
Other Names	CXCR3; GPR9; C-X-C chemokine receptor type 3; CXCR-3; CKR-L2; G protein-coupled receptor 9; Interferon-inducible protein 10 receptor; IP-10 receptor; CD183
Accession No.	Swiss-Prot: P49682NCBI Gene ID: 2833
Uniprot	P49682
GeneID	2833;
SDS-PAGE MW	37kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of K562 cells using CXCR-3 Polyclonal Antibody

Note: This product is for in vitro research use only