

Ob-R Polyclonal Antibody

Catalog No: #41688

Package Size: #41688-1 50ul #41688-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ob-R Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Ob-R Polyclonal Antibody detects endogenous levels of Ob-R protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Ob-R.
Target Name	Ob-R
Other Names	LEPR; DB; OBR; Leptin receptor; LEP-R; HuB219; OB receptor; OB-R; CD295
Accession No.	Swiss-Prot: P48357NCBI Gene ID: 3953
Uniprot	P48357
GeneID	3953;
SDS-PAGE MW	132kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa K562 MCF7 cells using Ob-R Polyclonal Antibody

Note: This product is for in vitro research use only