

# Mouse Transmembrane protein 163 (TMEM163) ELISA Kit



Catalog No: #EK6440

Orders: order@signalwayantibody.com

Package Size: #EK6440-1 48T #EK6440-2 96T

Support: tech@signalwayantibody.com

## Description

Product Name	Mouse Transmembrane protein 163 (TMEM163) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse (Mus musculus)
Other Names	DC29; DKFZp566N034; DKFZp666J217; FLJ30176; SV31;
Accession No.	Q8C996
Uniprot	Q8C996
GeneID	72160;
Storage	<p>The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.</p> <p>The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days, and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).</p>

## Application Details

Detect Range:Request Information

Sensitivity:Request Information

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

**Detection Method:**SandwichTest principle:This assay employs a two-site sandwich ELISA to quantitate TMEM163 in samples. An antibody specific for TMEM163 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyTMEM163 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TMEM163 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TMEM163 bound in the initial step. The color development is stopped and the intensity of the color is measured.**Product Overview:**SV31, a novel synaptic vesicle membrane protein of 31 kDa with six putative transmembrane helices that, according to its membrane topology and phylogenetic relation, may function as a vesicular transporter. The vesicular allocation is demonstrated by subcellular fractionation, heterologous expression, immunocytochemical analysis of brain sections and immunoelectron microscopy.

The protein is expressed in select brain regions and contained in subpopulations of nerve terminals that immunostain for the vesicular glutamate transporter 1 and the vesicular GABA transporter VGaT (vesicular amino acid transporter) and may attribute specific and as yet undiscovered functions to subsets of glutamatergic and GABAergic synapses.

---

Note: This product is for in vitro research use only