

Mouse Transmembrane 4 L6 family member 20 (TM4SF20) ELISA Kit

Catalog No: #EK6477

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #EK6477-1 48T #EK6477-2 96T

Description

Product Name	Mouse Transmembrane 4 L6 family member 20 (TM4SF20) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse (<i>Mus musculus</i>)
Other Names	FLJ22800; PRO994; TCCE518;
Accession No.	Q9CQY8
Uniprot	Q9CQY8
GeneID	66261;
Storage	<p>The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.</p> <p>The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days, and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).</p>

Application Details

Detect Range:Request Information

Sensitivity:Request Information

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:SandwichTest principle:This assay employs a two-site sandwich ELISA to quantitate TM4SF20 in samples. An antibody specific for TM4SF20 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyTM4SF20 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TM4SF20 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TM4SF20 bound in the initial step. The color development is stopped and the intensity of the color is measured.**Product Overview:**TM4SF20, Belongs to the L6 tetraspanin family, Members of the L6 family of membrane proteins, a branch of the tetraspanin superfamily, are overexpressed in tumor cells from many types of cancers. However, direct evidence of their oncogenic activity has not been previously shown.

this gene is expressed at low level, only 16.6% of the average gene in this release.The gene contains 4 different gt-ag introns. Transcription produces one mRNA. There are 3 validated alternative polyadenylation sites. The spliced mRNA putatively encodes a good protein, containing L6 membrane domain, some transmembrane domains, apparently vertebrate specific Proteins are expected to localize in various compartments.

Note: This product is for in vitro research use only