

MACC1 Polyclonal Antibody

Catalog No: #41698

Package Size: #41698-1 50ul #41698-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	MACC1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms
Specificity	MACC1 Polyclonal Antibody detects endogenous levels of MACC1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human MACC1.
Target Name	MACC1
Other Names	MACC1; Metastasis-associated in colon cancer protein 1; SH3 domain-containing protein 7a5
Accession No.	Swiss-Prot: Q6ZN28NCBI Gene ID: 346389
Uniprot	Q6ZN28
GeneID	346389;
SDS-PAGE MW	97kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of L929 cells using MACC1 Polyclonal Antibody

Note: This product is for in vitro research use only