

Human Transmembrane 4 L6 family member 19 (TM4SF19) ELISA Kit

Catalog No: #EK6480

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #EK6480-1 48T #EK6480-2 96T

Description

Product Name	Human Transmembrane 4 L6 family member 19 (TM4SF19) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human (Homo sapiens)
Other Names	OCTM4; OTTHUMP00000208460 tetraspan membrane protein OCTM4
Accession No.	Q96DZ7
Uniprot	Q96DZ7
GeneID	116211;
Storage	<p>The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.</p> <p>The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days, and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).</p>

Application Details

Detect Range:Request Information

Sensitivity:Request Information

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:SandwichTest principle:This assay employs a two-site sandwich ELISA to quantitate TM4SF19 in samples. An antibody specific for TM4SF19 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyTM4SF19 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TM4SF19 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TM4SF19 bound in the initial step. The color development is stopped and the intensity of the color is measured.**Product Overview:**TM4SF19 Belongs to the L6 tetraspanin family.Members of the L6 family of membrane proteins, a branch of the tetraspanin superfamily, are overexpressed in tumor cells from many types of cancers. However, direct evidence of their oncogenic activity has not been previously shown.

this gene is expressed at low level, only 16.6% of the average gene in this release.The gene contains 4 different gt-ag introns. Transcription produces one mRNA. There are 3 validated alternative polyadenylation sites. The spliced mRNA putatively encodes a good protein, containing L6 membrane domain, some transmembrane domains, apparently vertebrate specific Proteins are expected to localize in various compartments.

Note: This product is for in vitro research use only