

## Osteoglycin Polyclonal Antibody

Catalog No: #41700

Package Size: #41700-1 50ul #41700-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Osteoglycin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms
Specificity	Osteoglycin Polyclonal Antibody detects endogenous levels of Osteoglycin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Osteoglycin.
Target Name	Osteoglycin
Other Names	OGN; OIF; SLRR3A; Mimecan; Osteoglycin; Osteoinductive factor; OIF
Accession No.	Swiss-Prot: P20774NCBI Gene ID: 4969
Uniprot	P20774
GeneID	4969;
SDS-PAGE MW	35kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of L929 cells using Osteoglycin Polyclonal Antibody

---

Note: This product is for in vitro research use only