

VEGF-C Polyclonal Antibody

Catalog No: #41702

Package Size: #41702-1 50ul #41702-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VEGF-C Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	VEGF-C Polyclonal Antibody detects endogenous levels of VEGF-C protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human VEGF-C.
Target Name	VEGF-C
Other Names	VEGFC; Vascular endothelial growth factor C; VEGF-C; Flt4 ligand; Flt4-L; Vascular endothelial growth factor-related protein; VRP
Accession No.	Swiss-Prot: P49767NCBI Gene ID: 7424
Uniprot	P49767
GeneID	7424;
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of L929 cells using VEGF-C Polyclonal Antibody

Note: This product is for in vitro research use only