

IL-29 Polyclonal Antibody

Catalog No: #41703

Package Size: #41703-1 50ul #41703-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-29 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	IL-29 Polyclonal Antibody detects endogenous levels of IL-29 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human IL-29.
Target Name	IL-29
Other Names	IL29; IFNL1; ZCYTO21; Interleukin-29; IL-29; Cytokine Zcyto21; Interferon lambda-1; IFN-lambda-1
Accession No.	Swiss-Prot: Q8IU54NCBI Gene ID: 282618
Uniprot	Q8IU54
GeneID	282618;
SDS-PAGE MW	21kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of L929 cells using IL-29 Polyclonal Antibody

Note: This product is for in vitro research use only