

PDGF-D Polyclonal Antibody

Catalog No: #41704

Package Size: #41704-1 50ul #41704-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PDGF-D Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	PDGF-D Polyclonal Antibody detects endogenous levels of PDGF-D protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human PDGF-D.
Target Name	PDGF-D
Other Names	PDGFD; IEGF; SCDGFB; MSTP036; Platelet-derived growth factor D; PDGF-D; Iris-expressed growth factor; Spinal cord-derived growth factor B; SCDGF-B
Accession No.	Swiss-Prot: Q9GZP0NCBI Gene ID: 80310
Uniprot	Q9GZP0
GeneID	80310;
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of L929 cells using PDGF-D Polyclonal Antibody

Note: This product is for in vitro research use only