

Inhibin  $\beta$ -E Polyclonal Antibody

Catalog No: #41705

Package Size: #41705-1 50ul #41705-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Inhibin $\beta$ -E Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Inhibin $\beta$ -E Polyclonal Antibody detects endogenous levels of Inhibin $\beta$ -E protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Inhibin $\beta$ -E.
Target Name	Inhibin $\beta$ -E
Other Names	INHBE; Inhibin beta E chain; Activin beta-E chain
Accession No.	Swiss-Prot: P58166NCBI Gene ID: 83729
Uniprot	P58166
GeneID	83729;
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.

## Images

Western Blot analysis of L929 cells using Inhibin  $\beta$ -E Polyclonal Antibody

---

Note: This product is for in vitro research use only